WEEKLY MRI EQUIPMENT QUALITY CONTROL FOR SMALL PHANTOM

MR Facility Name: MR Scanner Identifier:

MR Facility Name: MR Scanner Identifier:													
1	2	3	4	5	6	7	8	9	10	11	12		
	Setup & Table Position		Center Freq	TX Gain or Attenua- tion	Measurements (Axial Slice #5 Diameter)			High-Contrast Spatial Resolution (Slice 1-Highest Resolved)		Detectability (Slice #)	Artifact Evaluation		
	Accuracy		(Hz)										
	Accuracy	Console		(dB)	H/F Sagittal Localizer	A/P	R/L	Upper	Lower	# of	Any	Tested	
Date	OK?	OK?			Length (mm)	(mm)	(mm)	Left	Right	Spokes	present?	Ву	Notes
Action limits:	±5 mm	Yes/No			100 ± 2 mm	100 ± 2 mm	100 ± 2mm	≤ 0.8 mm	≤ 0.8 mm		Yes/No		
	<u> </u>												
								_					

Reviewed by:	Date of Review:
· ————————————————————————————————————	