

Stereotactic Breast Biopsy Phantom Control Chart

Site: _____ Room: _____ Year: _____

Image Receptor: Film/Screen _____ Digital _____

Phantom Used: ACR Mammo Accreditation Phantom Digital "Mini" Phantom Other _____

Technique Factors: Target/Filter: _____ kVp Setting: _____

Phototimed Manual
Density Control Setting: _____ mA Setting: _____

Exposure Time (or mAs): _____

Phantom Control Chart – Phantom Scores

Month																								
Date																								
Performed by																								

Fibers	+2.0																								
	+1.0																								
	-1.0																								
	-2.0																								

Specks	+2.0																								
	+1.0																								
	-1.0																								
	-2.0																								

Masses	+2.0																								
	+1.0																								
	-1.0																								
	-2.0																								

Date:

Action:

Figure 5A: Phantom Control Chart.