

Stereotactic Breast Biopsy Hardcopy Output Quality Control

Site: _____ Room: _____ Year: _____

Hardcopy Display Unit: _____ Image(s) Used for Hardcopy Output Check: _____

Note: Circle Values Established as Control Values

Month		Jan	Feb	Mar	Apr	May	Jun	Jul	Aug	Sept	Oct	Nov	Dec
Filming Window Level													
Filming Window Width													
Hardcopy & Monitor Demonstrate Comparable Gray Scales (Yes/No)													
F I L M O D	Location #1 _____												
	Location #2 _____												
	Location #3 _____												
	Location #4 _____												
All measured Film ODs are within ± 0.2 of Control Levels (Yes/No)													
Check Performed By (initials)													
Physician Review *													
Medical Physicist Review *													

Date:

Action:

Figure 6. Stereotactic breast biopsy hardcopy output quality control form.