

**STEREOTACTIC BREAST BIOPSY
SCREEN-FILM REPEAT ANALYSIS SUMMARY**

Site: _____

Date From: _____ To: _____

Cause	Number of Films	Percentage Of Repeats
1. Positioning		
2. Patient Motion		
3. Light Films		
4. Dark Films		
5. Black Film		
6. Static		
7. Fog		
8. Incorrect Patient I.D., or Double Exposure		
9. Mechanical		
10. Artifacts		
11. Other		

$$\text{Repeat Rate} = \frac{\text{Total Repeated SBB}}{\text{Total SBB Films Used}} \times 100\%$$

Total Repeated SBB Films 1-11	
Total SBB Films Used	
Repeat Rate	

Performed by: _____ Date: _____

Action for Quality Improvement: _____

Physician Review _____
Medical Physicist Review _____

Date: _____
Date: _____

Figure 8. Stereotactic breast biopsy screen-film repeat analysis data collection