

**STEREOTACTIC BREAST BIOPSY
DIGITAL REPEAT ANALYSIS WORKSHEET**

(Enter any repeated exposures that required the patient to have additional dose beyond that of the normal exam.)

Site: _____

Stereotactic Breast Biopsy System: _____

Date From: _____ To: _____

| Causes | # of Times | Date | Technologist |
|--------|------------|------|--------------|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

Causes:

- | | |
|---|----------------------------|
| 1. Positioning | 6. X-Ray Equipment Failure |
| 2. Patient Motion | 7. Software Failure |
| 3. Detector Under exposure (Excessively noisy Images) | 8. Blank Image |
| 4. Detector Overexposure (saturation) | 9. Other |
| 5. Incorrect Patient ID | |

Figure 9. Stereotactic breast biopsy digital repeat analysis data collection worksheet